

Status: This is the original version (as it was originally made). This item of legislation is currently only available in its original format.

SCHEDULE 2

METHODS OF ANALYSIS

4b.

DETERMINATION OF TOTAL NITROGEN IN CALCIUM CYANAMIDE — IN THE PRESENCE OF NITRATE

Standard nitrate-ammoniacal solution:

4.13 Weigh out 2,500 g of potassium nitrate and 10.160 g of ammonium sulphate and place them in a 250 ml graduated flask. Dissolve in water and make up to 250 ml. 1 ml of this solution contains 0.010 g of nitrogen.