

SCHEDULE 2
TESTING METHODS

PART I

METHOD FOR THE ISOLATION OF *clostridium perfringens*

Inoculations

4. Once the agar has set, each agar plate shall be overlaid with a further 10 ml EY-free TSC agar at a temperature of $46^{\circ}\text{C} \pm 1^{\circ}\text{C}$. Once the overlay has set and with the plate lids uppermost the plates shall be incubated anaerobically at $37^{\circ}\text{C} \pm 1^{\circ}\text{C}$ for 20 hours \pm 2 hours.